Frederick National Laboratory for Cancer Research

sponsored by the National Cancer Institute



Annotation Visualization and Impact Analysis AVIA

Hue Vuong Reardon Advanced Biomedical Computational Science April 23, 2020


DEPARTMENT OF HEALTH AND HUMAN SERVICES • National Institutes of Health • National Cancer Institute

Overview

- Background
 - Sequencing and variants
 - Variant annotations
 - Impact analysis
- Demo
 - Submitting to AVIA and retrieving results
 - Single sample example
 - Multi sample example
 - Registered users
 - Project management
 - Cohort annotations

NGS Workflow





What is Annotation and Impact Analysis?




- Annotation: To process of identifying locations of genes and all of the coding regions in a genome.
 - Identifying other associated data for a variant at a given position
- Impact Analysis: To determine the unexpected, negative effects of a specific nucleotide change, or a set of nucleotide changes to the DNA of individuals
 - How do the changes affect individuals or populations?
 - What are the functions of the genes affected?

How do we sift through the hundreds/thousands of variants to find those of interest?







Can also be intergenic (between genes) or in a non-coding gene (ncRNA)


Gene Regulation



Methylation Promoters Tf Binding Sites

Splicing miRNA


PTM
Protein Binding
Binding Sites
Domains, etc



http://pubs.rsc.org/services/images/RSCpubs.ePlatform.Service.FreeContent.ImageService.svc/ImageService/Articleimage/2015/MB/c5mb00310e/c5mb00310e-f1_hi-res.gif







HMDD snoRNA miRNA


TargetScan SomamiR microPIR

Source: Curtan, A and Phillip Sharp. The Role of miRNAs in Regulating Gene Expression Networks. J. of Mol. Biol. (2013) 425(19):3582-3600.



Annotation and Functional Prediction

Now let's take a quick look at some ways of predicting and visualising the effect of variation on protein structure and function.





Annotation and Impact Analysis

- Annotation: Identifying other associated data at a variant's genomic location
 - Presence of gene or regulatory regions
 - Uniqueness and or repeat regions
 - Presence in other samples or studies
- Impact Analysis: Assessing the impact of that change
 - gene/protein/pathway
 - Pathogenicity predictions


How do we prioritize the hundreds/thousands of variants?

Annotations

- Gene RefSeq, (Ensembl)
- Regulatory regions TargetScan, HMDD,
- Population databases dbSNP, gnomAD, 1000 genomes
- Disease associated variants COSMIC, ClinVar, TCGA
- Genomic Features Genomicsuperdups, nonb, ENCODE
- Protein Features Prosite_domain, dbptm
- Protein scoring algorithms SIFT, polyphen, CADD
- 88 annotations in current version
- Regular updates through automated downloads

NGS Workflow






AVIA





All variants detected in sample(s)





AVIA Allele Frequency Designations

- We adopted these terms from Marth's Lab gene.iobio. The bins are as follows and annotations display the group with the highest allele frequency:
 - common AF >5%
 - uncommon AF 1-5%
 - rare AF = <1%
 - superrare <0.1%
 - uberrare AF <0.01%






Reference Genomes in AVIA

- Human
 - UCSC hg19 NCBI GRCh37 (current)
 - UCSC hg38 NCBI GRCh38
- Mouse
 - UCSC mm10 GRCm38





Variant Call Format (VCF) is the preferred format



Format and Sample go hand in hand

Variant Types

- Single base-pair substitution
 - Single nucleotide polymorphisms (SNPs)
- Multiple nucleotide substitution
 - Substitutions where length > 1
- Insertion or deletion, also known as 'indel'
 - Insertion or deletion of a DNA sequence
 - 2 to 100's of base-pairs in length For AVIA, limited to small indels < 50
- Structural variation
 - larger DNA sequence
 - copy number variation
 - chromosomal rearrangement events



Indel Representation




Fig. 1. Example of VCF entries representing the same variant. Left panel aligns each allele to the reference genome, and the right panel represents the variant in VCF. (**A**) is not left-aligned (**B**) is neither left-aligned nor parsimonious, (**C**) is not parsimonious and (**D**) is normalized



AVIA Indel Normalization

All indels are normalized using U. Michigan's VT package



- Annotations against normalized indels
- Indel alias table
 - Maintain all aliases













Annotation, Visualization, and Impact Analysis								
Home And	alysis∙	Examples+	About-	Projects+			Hello huetogo	Sign out
AVIAv3 A	Annot	ation Dat	What's FAQs AVIA Date	new dabase Sources	n(hg19)		Human GRCh38/hg38	Mouse (mm10)
Category 👢	Database	Name		↑ Version	↓↑ Description ↓	Last Updated	Search: Citation	11
Alternative Splicing	ALT_SP	LICE		NA	Ensembl Splice events	21-FEB- 17	Koscielny G, Le Texier Gopalakrishnan C, Ku Riethoven JJ, Nardone Fallsehr C, Hofmann C Harrington E, Boue S, M, Lopez F, Ritchie W, Ara T, Pospisil H, Herrr Reich J, Guigo R, Bork MK, Vilo J, Hide W, Ap Thanaraj TA, Gauther	manduri V, e F, Stanley E,), Kull M, Eyras E, Plass Moucadel V, nann A, G P, Doeberitz weiler R,
Disease Related	CANDL			20161222	Cancer Driver Log (CanDL): Catalog of Potentially	21-FEB-		

https://avia-abcc.ncifcrf.gov



Impact Assessment

- Variant overview, analytics
- Gene gene.iobio
- Protein ProtVista, MolArt
- Gene Functional clustering DAVID
- Pathway PathView
- Tissue SAMM

- Literature references
- Comparisons between and within annotations, samples



AVIA Demo Overview

- Basic Navigation
- Submit variant list to AVIA
- Data Retrieval
- Walk through of visualization and data
- Advanced Features
- Registration and Additional Tools
 - Custom Annotations
 - Project Management
 - Data Sharing
 - Saving and sharing dashboards
 - Building cohorts
 - Reannotating



Basic Navigation

Navigation








Annotation, Visualization, and Impact Analysis

Home Analysis Examples About Sign in

AVIA v3.0 Features




URL: https://avia-abcc.ncifcrf.gov






Use any tool or retrieve results



View a sample analysis










Email: NCIFrederickAVIA@mail.nih.gov

Web form to ask questions or request services





Information about AVIA





Accepted Input Formats



Understanding VCF format

- VCF variant call format
- standard file format for storing variation data
- used by large scale variant mapping projects
- the standard output of variant calling software
- can be compressed and indexed

VCF is a preferred format because it is unambiguous, scalable and flexible, allowing extra information to be added to the info field. Many millions of variants can be stored in a single VCF file.




Header in VCF files

```
##fileformat=VCFv4.3
##fileDate=20090805
##source=myImputationProgramV3.1
##reference=file:///seq/references/1000GenomesPilot-NCBI36.fasta
##contig=<ID=20,length=62435964,assembly=B36,md5=f126cdf8a6e0c7f379d618ff66beb2da,species="Homo sapiens",taxonomy=x>
##phasing=partial
##INFO=<ID=NS, Number=1, Type=Integer, Description="Number of Samples With Data">
##INFO=<ID=DP, Number=1, Type=Integer, Description="Total Depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele Frequency">
##INFO=<ID=AA, Number=1, Type=String, Description="Ancestral Allele">
##INFO=<ID=DB, Number=0, Type=Flag, Description="dbSNP membership, buard 129">
##INFO=<ID=H2, Number=0, Type=Flag, Description="HapMap2 membership">
##FILTER=<ID=q10,Description="Quality below 10">
##FILTER=<ID=s50,Description="Less than 50% of samples have data">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=GO, Number=1, Type=Integer, Description="Genotype Quality">
##FORMAT=<ID=DP, Number=1, Type=Integer, Description="Read Depth">
##FORMAT=<ID=HQ, Number=2, Type=Integer, Description="Haplotype Quality">
#CHROM POS
                           REF
                                 ALT
                                         OUAL FILTER
                                                       INFO
                                                                                          FORMAT
                                                                                                       NA00001
                                                                                                                        NA00002
                                                                                                                                         NA00003
       14370
                rs6054257 G
                                                       NS=3;DP=14;AF=0.5;DB;H2
20
                                               PASS
                                                                                          GT:GQ:DP:HQ
                                                                                                       0 0:48:1:51,51 1 0:48:8:51,51
1/1.43.5....
       17330
                                               a10
                                                       NS=3:DP=11:AF=0.017
                                                                                          GT:GQ:DP:HQ 0|0:49:3:58,50 0|1:3:5:65,3
                                                                                                                                         0/0:41:3
20
       1110696 rs6040355 A
                                         67
                                                       NS=2;DP=10;AF=0.333,0.667;AA=T;DB GT:GQ:DP:HQ 1|2:21:6:23,27 2|1:2:0:18,2
                                                                                                                                         2/2:35:4
20
                                               PASS
                                                                                                                                         0/0:61:2
       1230237
                                               PASS
                                                       NS=3;DP=13;AA=T
                                                                                                       0|0:54:7:56,60 0|0:48:4:51,51
20
                                                                                          GT:GQ:DP:HQ
20
       1234567
                microsat1 GTC
                                 G,GTCT 50
                                               PASS
                                                       NS=3;DP=9;AA=G
                                                                                          GT:GO:DP
                                                                                                       0/1:35:4
                                                                                                                        0/2:17:2
                                                                                                                                         1/1:40:3
```



VCF (ctd)

VCF files are tab-delimited text files



position and alt allele

Contains info about variant e.g. counts, allele frequencies (AF), depth, etc Sample Information

1 – many samples in a single VCF file

Format and Sample go hand in hand



Sample Columns in VCF file

```
##FILTER=<ID=SS0,Description= Less than S0% of Samples have data >
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
##FORMAT=<ID=GQ,Number=1,Type=Integer,Description="Genotype Quality">
##FORMAT=<ID=DP,Number=1,Type=Integer,Description="Read Depth">
##FORMAT=<ID=HQ,Number=2,Type=Integer,Description="Haplotype Quality">
#CHROM POS ID REF ALT QUAL FILTER INFO
```

```
FORMAT
                             NA00002
             NA00001
                                               NA00003
GT:GQ:DP:HQ
             0 0:48:1:51,51 1 0:48:8:51,51
             0 0:49:3:58,50
                             0 1:3:5:65,3
GT:GQ:DP:HQ
                                               0/0:41:3
GT:GQ:DP:HQ
             1 2:21:6:23,27
                             2 1:2:0:18,2
                                               2/2:35:4
                             0|0:48:4:51,51
             0 0:54:7:56,60
GT:GQ:DP:HQ/
                                               0/0:61:2
GT:GQ:DP
             0/1:35:4
                             0/2:17:2
                                               1/1:40:3
```

/ : genotype unphased

: genotype phased



BED – like (ANNOVAR format)

The first five space or tab delimited fields are Chromosome ("chr" prefix is optional), Start, End, Reference Allelel, Alternative Allele. The rest of the columns are completely optional.


```
cat ex1.avinput
                                        comments: rs15842, a SNP in 5' UTR of ISG15
        948921 948921 T
                                        comments: rs149123833, a SNP in 3' UTR of ATAD3C
        1404001 1404001 G
                                        comments: rs1287637, a splice site variant in NPHP4
        5935162 5935162 A
        162736463
                        162736463
                                                        comments: rs1000050, a SNP in Illumina SNP arrays
                                                        comments: rs6576700 or SNP A-1780419, a SNP in Affymetrix SNP arrays
        84875173
                        84875173
        13211293
                        13211294
                                                        comments: rs59770105, a 2-bp deletion
        11403596
                        11403596
                                                AT
                                                        comments: rs35561142, a 2-bp insertion
                                                        comments: rs10552169, a block substitution
        105492231
                                                ATAAA
                        105492231
        67705958
                        67705958
                                                        comments: rs11209026 (R3810), a SNP in TL23R associated with Crohn's disease
                                                         comments: rs2241880 (T300A), a SNP in the ATG16L1 associated with Crohn's disease
        234183368
                        234183368
        50745926
                        50745926
                                                        comments: rs2066844 (R702W), a non-synonymous SNP in NOD2
        50756540
                        50756540
                                                        comments: rs2066845 (G908R), a non-synonymous SNP in NOD2
        50763778
                        50763778
                                                        comments: rs2066847 (c.3016 3017insC), a frameshift SNP in NOD2
13
                                                        comments: rs1801002 (del35G), a frameshift mutation in GJB2, associated with hearing loss
        20763686
                        20763686
                                                        comments: a 342kb deletion encompassing GJB6, associated with hearing loss
        20797176
                        21105944
        8887543 8887543 A
                                        comments: a mutation that abolishes stop codon
        8887539 8887539 A
                                        comments: a mutation that results in premature stop codon
                                GATT
                                        comments: a mutation that creates a stop codon 2 amino acids downstream
        8887536 8887537 AG
                                GGAA
                                        comments: a mutation that results in insertion of a new amino acid
        8887540 8887540 G
        1295288 1295288 G
                                        comments: a variant upstream of transcriptional start site
        95602958
                        95602958
                                                        comments: a variant that affects splicing of UTR regions
```

Indels cannot be normalized using this input format. VCF format is the preferred method.

https://doc-openbio.readthedocs.io/projects/annovar/en/latest/user-guide/input/







Refer to ANNOVAR about normalization and variant highjacking: https://doc-openbio.readthedocs.io/projects/annovar/en/latest/articles/VCF/

Sample Mining for Single Sample (ss VCF, BED, or dbSNP ids)



- AVIA does not mine BED files for sample information
- When no sample information is presented, a dummy sample id is used instead



Submit to AVIA








Accepted Input Types

- VCF (single and multisample)
- Bed-like (Chr start stop ref alt)
- dbSNP identifiers

AVIA Annotation and Visualization Request









Your results for vizssbed can be found at <u>AVIA</u>. You can get your results one of two ways: 1) If you are not a registered user, click on the 'Retrieve Results' under the 'Analysis' menu on the link above. Then enter your AVIA id <u>vizssbed</u> and email address. 2) If you are a registered user, please log into your AVIA Account by selecting "Sign in" and follow directions for CILogon. After authentication, select your AVIA id from your personal dashboard. Thank you for using AVIA. If the link does not work, please use https://avia-abcc.ncifcrf.gov.



1)

3)

avia-abcc.ncifcrf.gov





Annotation, Visualization, and Impact Analysis

2) Analysis

Examples- About-




AVIA v3.0 Features

AVIAv3 Annotation Databases for Human

Analysis - Examples
Analyze Variant

Retrieve Results

Customized Tools



Results Navigation Landing Page







Demo



Questions?

Contact Us:

NCI-FrederickAVIA@mail.nih.gov

Or https://avia-abcc.ncifcrf.gov